

Product Catalog

Rapid Test Kits & Uncut Sheets (Version 1.06)



Biogate Laboratories Ltd.

110-4238 Lozells Avenue,
Burnaby, BC, Canada,
V5A 0C4

Phone: 604-322-2955

Fax: 604-322-2995

E-mail: info@biogatelab.com

Biogate Laboratories Ltd. BC, Canada

www.biogatelab.com

Biogate Laboratories Ltd. was founded in 2004 by a team of talented scientists in Vancouver, Canada. We are a research – orientated company and dedicated to developing, producing and marketing primary and secondary antibodies, purified antigens and proprietary blocking reagents.

After years development, our mature technology enables us to offer a comprehensive product line in the rapid test field. Currently, our colloidal gold rapid test product line has expanded to more than twenty major products, majorly covering all type of the disease in the global market.

We believe client's success is our success. With a commitment providing high quality products and service to our partners, we will continuously create more mutual benefits for both of us in the future.

Rapid Diagnostic Experts From Canada Biogate Laboratories Ltd.

Add: 110-4238 Lozells Avenue, Burnaby,
Canada, V5A 0C4
TEL: 1-604-322-2955
FAX: 1-604-322-2995
EMAIL: info@biogatelab.com

Table of Contents

Fertility Hormone	3
TORCH	4
Hepatitis Disease	5
Sexually Transmitted Diseases	6
Tropic Disease	7
Gastrointestinal Infection	8
Cardiac Marker	8
Drug of Abuse Rapid Test	9
Cancer Marker and Other Infection Disease	10



Cancer Marker & Other Disease Rapid Test Kits (Colloidal Gold)*

Catalog	Product	Format**	Specimen	Kit Size
RCA-010C	PSA	Cassette	S/P	20 test
RB2-010C	β 2-Microglobulin	Cassette	Urine	20 test
RB2-110C	β 2-Microglobulin	Cassette	Serum	20 test
RCA-020C	Fecal Occult Blood	Cassette	Stool	20 test
RCP-017C	C. Pneumoniae IgM	Cassette	S/P	20 test
RCP-018C	C. Pneumoniae IgG	Cassette	S/P	20 test
RTB-020C	Tuberculosis Ab	Cassette	S/P	20 test
RMP-017C	M.Pneumonia IgM	Cassette	S/P	20 test
RMP-018C	M.Pneumonia IgG	Cassette	S/P	20 test
RCT-010C	C. Trachomatis	Cassette	S/P	20 test

* corresponding uncut sheets for different application are available

** cassette format is available



HCG Rapid Test Kits (Colloidal Gold)

Catalog	Product	Format	Specimen	Sensitivity	Kit Size
RFH-010S	HCG	Strip	Urine	25 mIU/mL	50 test
RFH-010C	HCG	Cassette	Urine	25 mIU/mL	25 test
RFH-010M	HCG	Midstream	Urine	25 mIU/mL	20 test
RFH-110S	HCG Combo	Strip	Urine/Serum	25 mIU/mL	50 test
RFH-110C	HCG Combo	Cassette	Urine/Serum	25 mIU/mL	25 test
RFH-110M	HCG Combo	Midstream	Urine/Serum	25 mIU/mL	20 test

LH Rapid Test Kits (Colloidal Gold)

Catalog	Product	Format	Specimen	Sensitivity	Kit Size
RFH-020S	LH	Strip	Urine	25 mIU/mL	50 test
RFH-020C	LH	Cassette	Urine	25 mIU/mL	25 test
RFH-020M	LH	Midstream	Urine	25 mIU/mL	20 test

HCG/LH Uncut Sheet Product

Catalog	Product	Application*	Specimen	Sensitivity	Size**
UFH-010S	HCG	Strip	Urine	25 mIU/mL	300 mm x W
UFH-010C	HCG	Cassette	Urine	25 mIU/mL	300 mm x W
UFH-110S	HCG	Strip	Urine/Serum	-	300 mm x W
UFH-110C	HCG	Cassette	Urine /Serum	-	300 mm x W
UFH-020S	LH	Strip	Urine	30 mIU/mL	300 mm x W
UFH-020C	LH	Cassette	Urine	30 mIU/mL	300 mm x W

* midstream format is available (UFH-010M/UFH-110M)

** width (W) of all kinds of uncut sheet is customized

TORCH Rapid Test Kits (Colloidal Gold)


Catalog	Product	Format	Specimen *	Kit Size
RTO-016C	Toxoplasma IgG/IgM	Cassette	S/P	25 test
RTO-026C	Rubella Virus IgG/IgM	Cassette	S/P	25 test
RTO-036C	CMV IgG/IgM	Cassette	S/P	25 test
RTO-046C	HSV-2 IgG/IgM	Cassette	S/P	25 test

TORCH Uncut Sheet Products

Catalog	Product	Application*	Specimen *	Size**
UTO-016C	Toxoplasma IgG/IgM	Cassette	S/P	300 mm x W
UTO-026C	Rubella Virus IgG/IgM	Cassette	S/P	300 mm x W
UTO-036C	CMV IgG/IgM	Cassette	S/P	300 mm x W
UTO-046C	HSV-2 IgG/IgM	Cassette	S/P	300 mm x W

* whole blood test could be conducted with assistant of assay buffer; but it may not work for every patient due to blood coagulation .

** width (W) of all kinds of uncut sheet is customized



Drug of Abuse Rapid Test Kits (Colloidal Gold)*

Catalog	Product	Format**	Specimen	Kit Size
RDA-010S	Amphetamines	Strip	Urine	40 test
RDA-020S	Barbital	Strip	Urine	40 test
RDA-030S	Benzodiazepines	Strip	Urine	40 test
RDA-040S	Cocaine	Strip	Urine	40 test
RDA-050S	Ecstasy	Strip	Urine	40 test
RDA-060S	Methamphetamines	Strip	Urine	40 test
RDA-070S	Methadone	Strip	Urine	40 test
RDA-080S	Morphine	Strip	Urine	40 test
RDA-090S	Opiates/Heroin/Morphine	Strip	Urine	40 test
RDA-100S	Marijuana	Strip	Urine	40 test
RDA-110S	Tricyclic Antidepressant	Strip	Urine	40 test
RDA-120S	Buprenorphine	Strip	Urine	40 test
RDA-130S	EDDP	Strip	Urine	40 test
RDA-140S	Ketamine	Strip	Urine	40 test
RDA-150S	Oxycodone	Strip	Urine	40 test
RDA-160S	Propoxyphene	Strip	Urine	40 test
RDA-170S	Phencyclidine	Strip	Urine	40 test

* corresponding uncut sheets for different application are available , please ask for more details

** cassette format is available

Gastrointestinal Infection Disease

Gastrointestinal Infection Rapid Test Kits/ Uncut Sheets (Colloidal Gold)

Catalog	Product	Format	Specimen *	Kit Size
RGI-016C	H.Pylori IgG/IgM	Cassette	S/P	20 test
RGI-011C	H.Pylori Ag	Cassette	Stool	20 test
Catalog	Product	Application	Specimen *	Size
UGI-011C	H.Pylori Ag	Cassette	Stool	300 mm x W
UGI-016C	H.Pylori IgG/IgM	Cassette	S/P	300 mm x W

Cardiac Marker

Cardiac Marker Rapid Test Kits/ Uncut Sheets (Colloidal Gold)

Catalog	Product	Format	Specimen *	Kit Size
RCM-010C	Troponin I	Cassette	S/P	20 est
RCM-020C	CK-MB	Cassette	S/P	20 test
RCM-030C	Human Myoglobin	Cassette	S/P	20 test
Catalog	Product	Application	Specimen *	Size**
UCM-010C	Troponin I	Cassette	S/P	300 mm x W
UCM-020C	CK-MB	Cassette	S/P	300 mm x W
UCM-030C	Human Myoglobin	Cassette	S/P	300 mm x W

* whole blood test could be conducted with assistant of assay buffer; but it may not work for every patient due to blood coagulation

** width (W) of all kinds of uncut sheet is customized

Hepatitis Disease

Hepatitis Disease Rapid Test Kits (Colloidal Gold)

Catalog	Product	Format	Specimen *	Kit Size
RHE-016C	HAV IgG/IgM	Cassette	S/P	25 test
RHE-021S	HBsAg	Strip	S/P	50 test
RHE-021C	HBsAg	Cassette	S/P	25 test
RHE-031S	HBeAg	Strip	S/P	50 test
RHE-031C	HBeAg	Cassette	S/P	25 test
RHE-049C	HCV	Strip	S/P	50 test
RHE-049C	HCV	Cassette	S/P	25 test
RHE-056C	HEV IgG/IgM	Cassette	S/P	25 test

Hepatitis Disease Uncut Sheet Products

Catalog	Product	Application*	Specimen *	Size**
UHE-016C	HAV IgG/IgM	Cassette	S/P	300 mm x W
UHE-020S	HBsAg	Strip	S/P	300 mm x W
UHE-021C	HBsAg	Cassette	S/P	300 mm x W
UHE-031S	HBeAg	Strip	S/P	300 mm x W
UHE-031C	HBeAg	Cassette	S/P	300 mm x W
UHE-049S	HCV	Strip	S/P	300 mm x W
UHE-049C	HCV	Cassette	S/P	300 mm x W
UHE-056C	HEV IgG/IgM	Cassette	S/P	300 mm x W

* whole blood test could be conducted with assistant of assay buffer; but it may not work for every patient due to blood coagulation

** width (W) of all kinds of uncut sheet is customized

Sexually Transmitted Disease

Sexually Transmitted Diseases (STD) Rapid Test Kits (Colloidal Gold)

Catalog	Product	Format	Specimen *	Kit Size
RSD-019S	TP Syphilis Ab	Strip	S/P	40 test
RSD-019C	TP Syphilis Ab	Cassette	S/P	20 test
RSD-016C	Syphilis IgG/IgM	Cassette	S/P	20 test
RSD-121C	HIV 1/2 Tri-line	Cassette	S/P	20 test
RSD-021C	HIV 1/2 Duo-line	Cassette	S/P	20 test

Sexually Transmitted Diseases (STD) Uncut Sheets Products

Catalog	Product	Application	Specimen *	Size**
USD-019S	TP Syphilis Ab	Strip	S/P	300 mm x W
USD-019C	TP Syphilis Ab	Cassette	S/P	300 mm x W
USD-016C	Syphilis IgG/IgM	Cassette	S/P	300 mm x W
USD-121C	HIV 1/2 Tri-line	Cassette	S/P	300 mm x W
USD-021C	HIV 1/2 Duo-line	Cassette	S/P	300 mm x W

* whole blood test could be conducted with assistant of assay buffer; but it may not work for every patient due to blood coagulation

** width (W) of all kinds of uncut sheet is customized

Tropical Disease

Tropical Disease Rapid Test Kits (Colloidal Gold)

Catalog	Product	Format	Specimen *	Kit Size
RDN-011C	Dengue NS1	Cassette	WB/S/P	20 test
RDN-016C	Dengue IgG/IgM	Cassette	S/P	20 test
RDN-013C	Dengue NS1/IgG/IgM	Cassette	WB/S/P	15 test
RMA-019C	Malaria pv/pf Ab	Cassette	S/P	20 test
RMA-011C	Malaria pf/pan Ag	Cassette	S/P	20 test
RST-016C	S. Typhoid IgG/IgM	Cassette	WBS/P	20 test
RST-011C	S. Typhoid Ag	Cassette	Stool/S/P	20 test

Tropical Disease Uncut Sheets

Catalog	Product	Application	Specimen *	Size
UDN-011C	Dengue NS1	Cassette	WB/S/P	300 mm x W
UDN-016C	Dengue IgG/IgM	Cassette	S/P	300 mm x W
UMA-011C	Malaria pf/pan Ag	Cassette	S/P	300 mm x W
UMA-019C	Malaria pv/pf Ab	Cassette	S/P	300 mm x W
UST-011C	S. Typhoid Ag	Cassette	Stool/S/P	300 mm x W
UST-016C	S. Typhoid IgG/IgM	Cassette	WB/S/P	300 mm x W

* whole blood test could be conducted with assistant of assay buffer; but it may not work for every patient due to blood coagulation

** width (W) of all kinds of uncut sheet is customized